

## rhodotorulapepsin

Cat. No. EXWM-4264

Lot. No. (See product label)

## Introduction

**Description** From the imperfect yeast Rhodotorula glutinis. Somewhat similar enzymes have been isolated from the

imperfect yeast-like organism Cladosporium sp. and the imperfect fungus Paecilomyces varioti.

**Synonyms** Rhodotorula aspartic proteinase; Cladosporium acid protease; Cladosporium acid proteinase;

Paecilomyces proteinase; Cladosporium aspartic proteinase; Paecilomyces proteinase; Rhodotorula glutinis aspartic proteinase; Rhodotorula glutinis aspartic

II; Rhodotorula acid proteinase

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 3.4.23.26

*CAS No.* 37259-59-9

**Reaction** Specificity similar to that of pepsin A. Cleaves Z-Lys+Ala-Ala and activates trypsinogen

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

## Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C $\sim$ -80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1