

teichoic-acid-transporting ATPase

Cat. No. EXWM-4679

Lot. No. (See product label)

Introduction

Description Order Teichoic Acid-Transporting ATPase for bacterial cell wall studies. Creative Enzymes provides high-

quality enzymes for research.

Product Information

Form Liquid or lyophilized powder

EC Number EC 7.5.2.4 (Formerly EC 3.6.3.40)

Reaction ATP + H2O + teichoic acidin = ADP + phosphate + teichoic acidout

Notes This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1