

Halohydrin Dehalogenase (Crude Enzyme)

Cat. No. NATE-1844

Lot. No. (See product label)

Introduction

Description This enzyme belongs to the family of hydrolases, specifically those acting on halide bonds in carbon-halide compounds. This product with the indicated enzyme activity was briefly purified from engineered E. coli.

Applications environmental protection; synthesis; analysis

Synonyms 2-haloacid dehalogenase[ambiguous]; 2-haloacid halidohydrolase [ambiguous][ambiguous]; 2-haloalkanoic acid dehalogenase; 2-haloalkanoic acid halidohydrolase; 2-halocarboxylic acid dehalogenase II; DL-2-haloacid dehalogenase[ambiguous]; L-2-haloacid dehalogenase; L-DEX

Product Information

Source E. coli

Appearance Clear to translucent yellow solution

EC Number EC 3.8.1.2

CAS No. 37289-39-7

Activity Undetermined

Reaction (S)-2-haloacid + H₂O = (R)-2-hydroxyacid + halide

Notes Since this product needs to be freshly prepared, it will take about 2 weeks after you confirm the order. Each time of the freeze-thawing may cause partial inactivation. Therefore, it should be dispensed as required and stored at -20 °C or lower. With the preservation of the extension of time, the enzyme activity will decline to a certain extent, so the product should be used as soon as possible. This product may have turbidity or precipitation in the production and preservation process, it can be mixed after melting and will not affect the normal use. This product is limited to scientific research use, shall not be used for clinical diagnosis or treatment, shall not be used for food or medicine, shall not be stored in ordinary residential. For your safety and health, please wear an experimental suit and wear disposable gloves.

Usage and Packaging

Package 100ml

Storage and Shipping Information

Storage at -20 °C or lower, for at least 1 month.