

## Protease Inhibitor Cocktail II

*Cat. No.* CEI-0796

*Lot. No.* (See product label)

### *Introduction*

**Description** Inhibits aspartic, cysteine, serine, and metalloproteases, as well as aminopeptidases. Recommended for use with bacterial cell extracts. Reconstitute each vial with 1ml DMSO and 4ml H<sub>2</sub>O to obtain a 5ml stock solution containing 20 mM AEBSF.HCl, 1.7mM bestatin, 200uM E-64, 85mM EDTA, and 2mM pepstatin A.

### *Product Information*

**Appearance** Lyophilized solid

**Targets** Protease