

## Pefachrome Try

Cat. No. CSUB-0595

Lot. No. (See product label)

### Introduction

**Applications** Pefachrome substrates are used for chromogenic assays of various serine proteinases. Pefachrome Try is used for the detection of serine proteases such as trypsin.

**Synonyms** Bz-Val-Gly-Arg-4-nitroanilide acetate; N-Benzoyl-L-valyl-glycyl-L-arginine-4-nitroanilide acetate

### Product Information

**CAS No.** 161572-28-7

**Molecular Formula** C<sub>26</sub>H<sub>34</sub>N<sub>8</sub>O<sub>6</sub> · C<sub>2</sub>H<sub>4</sub>O<sub>2</sub>

**Molecular Weight** 614.65

**Substrates** Trypsin Substrates

### Storage and Shipping Information

**Storage** 2-8°C