

## tyrosine-ester sulfotransferase

Cat. No. EXWM-3403

Lot. No. (See product label)

### Introduction

**Description** Phenols and organic hydroxylamines can act as acceptors (cf. EC 2.8.2.1 aryl sulfotransferase).

**Synonyms** aryl sulfotransferase IV; L-tyrosine methyl ester sulfotransferase

### Product Information

**Form** Liquid or lyophilized powder

**EC Number** EC 2.8.2.9

**CAS No.** 9055-56-5

**Reaction** 3'-phosphoadenylyl sulfate + L-tyrosine methyl ester = adenosine 3',5'-bisphosphate + L-tyrosine methyl ester 4-sulfate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

### Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.