

## XTP/dITP diphosphatase

Cat. No. EXWM-4639

Lot. No. (See product label)

## Introduction

**Description** The enzymes from the bacterium Escherichia coli and the archaea Methanococcus jannaschii and

Archaeoglobus fulgidus are highly specific for XTP and dITP. The activity is dependent on divalent cations,

Mg2+ is preferred.

**Synonyms** hypoxanthine/xanthine dNTP pyrophosphatase; rdgB (gene name)

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 3.6.1.66

**Reaction** (1) XTP + H2O = XMP + diphosphate; (2) dITP + H2O = dIMP + diphosphate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

## Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1