

## elisabethatriene synthase

*Cat. No.* EXWM-5200

*Lot. No.* (See product label)

### **Introduction**

**Description** Requires Mg<sup>2+</sup> or less efficiently Mn<sup>2+</sup>. The enzyme is also able to use farnesyl diphosphate and geranyl diphosphate.

**Synonyms** elisabethatriene cyclase

### **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 4.2.3.41

**CAS No.** 334022-59-2

**Reaction** geranylgeranyl diphosphate = elisabethatriene + diphosphate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

### **Storage and Shipping Information**

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.