

## Enzyme stabilizer AES

*Cat. No.* AES-001

*Lot. No.* (See product label)

### **Introduction**

**Description** Broad-spectrum stabilizer for enzyme, antigen and antibody; it can be used in diagnostic reagents and enzymes diluent solutions, especially for various oxidases, glycerophosphate oxidases, and uricase, etc. Better effect when combined with HPD.

**Synonyms** Enzyme stabilizer AES; Enzyme stabilizer; AES

### **Product Information**

**Form** White powder

**Purity** ≥98.0%

### **Usage and Packaging**

**Package** 100g/Bottle

### **Storage and Shipping Information**

**Storage** Store at 2° C~8° C. Protect from moisture.