

Urea solution

Cat. No. CSUB-0608

Lot. No. (See product label)

Introduction

Applications Recommended for the detection of urease activity.

Synonyms Carbamide; Carbonyldiamide

Product Information

CAS No. 57-13-6

Molecular Weight 60.06

Concentration 40 % (w/v) in H₂O

Substrates Urease

Usage and Packaging

Preparation Instructions Warm the refrigerated 40% urea solution to room temperature. Aseptically add 5 mL to 95 mL of sterile molten Urea Agar Base (Christensen). Mix well and dispense in sterile tubes.

Storage and Shipping Information

Storage 2-8°C